










THUASNE

UniReliever™

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Properties/Mode of action

Adjustable 3-point leverage system with inelastic straps to provide alignment support and off-loading on the affected compartment.

This realignment helps reduce varus or valgus angles of the knee and helps distribute excessive load away from the damaged knee compartment.

Through correction of alignment, the brace will help reduce load, aiding in the reduction of compression and inflammation.

Reducing the load on the affected side would lead to the reduction of pain and cartilage degeneration.

The brace maintains proper fit on the leg through:

- inner silicon paddings (shells and straps),

- semi-rigid femoral and tibial shells allowing an intimate fit to the patient's unique morphology.

Easy to use thanks to:

- an offloading system with self-dosing adjustment,

- color coded clip buckles.

Extension can be adjusted to 0°, 5°, 10°, 15°, 20° and 30°.

Flexion can be adjusted to 0°, 30°, 45°, 60°, 75°, 90° and 110°.

Indications

Symptomatic unicompartimental femorotibial osteoarthritis (moderate to severe).

Knee off-loading for post-traumatic, post-operative or degenerative conditions. Joint instability/laxity.

Alternative to osteotomy or leg misalignment surgery.

Contraindications

Do not use the product if the diagnosis has not been confirmed.

Do not use in the event of known allergy to any of the components.

Do not apply the product in direct contact with broken skin or an open wound without an adequate dressing.

Do not use in the event of severe varicose veins (preventing the regular use of an offloading brace).

Do not use in the event of significant genu recurvatum.

Do not use in the event of significant genu varum or genu valgum.

History of venous or lymphatic disorders.

Do not use in the event of major venous thromboembolic history without thromboprophylaxis.

Precautions

Verify the product's integrity before every use.

Do not use the product if it is damaged.

It is recommended that a healthcare professional supervises the first application.

Strictly comply with your healthcare professional's prescription and recommendations for use.

This product is intended for the treatment of a given condition. Its duration of use is limited to this treatment only.

For hygiene and performance reasons, do not re-use the product for another patient.

It is recommended to adequately tighten the device to achieve support/immobilisation without restricting blood circulation.

In the event of discomfort, significant hindrance, pain, variation in limb volume, abnormal sensations or change in colour of the extremities, remove the device and consult a healthcare professional.

In the event of any modification in the product's performance, remove it and consult a healthcare professional.

Before any sports activity, check the compatibility of the use of this medical device with your healthcare professional.

Do not wear the product in a medical imaging machine.

Do not use the device in case of application of certain products on the skin (creams, ointments, oils, gels, patches,...).

Undesirable side-effects

This device can cause skin reactions (redness, itching, burns, blisters, etc.) or wounds of various degrees of severity. Any serious incidents related to the product should be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is resident.

Instructions for use/Application

It is recommended to wear the device directly on the skin, unless contraindicated. Please read these instructions with the patient when fitting the brace for the first time and make sure that he/she understands how to position the knee brace.

How to fit the knee brace:

• Step 1: Evaluate the patient

Assess the level of pain felt by the patient to measure the function and efficacy of the knee brace after it is fitted.

• Step 2: Positioning the brace on the leg

- Unclip and loosen all the straps of the knee brace, and fold the self-fastening straps back for easier application.

Ensure that the dial is well pulled so that there is no tension in the straps.

- Have the patient sit on the edge of a chair with their knee bent to approximately 30-60° of flexion.

- Place the knee brace on the bare leg (Fig. 1), the hinge on affected compartment side.

Align the hinge at knee center, corresponding with the upper third of your patella (Fig. 2).

- Close the lower calf strap buckle A, followed by the top thigh strap buckle B (Fig. 3).

Follow numbers and colors indications on the buckles and the shells.

This should secure the knee brace for positioning of the strap system, ensuring the hinge is positioned correctly.

Position the Cross Buttress (where the straps intersect) near the edge of patella (Fig. 4).

Use the marking on the Cross Buttress to ensure proper alignment with the patella (Fig. 5).

- Close front tibia strap buckle C, followed by the front thigh strap buckles D (Fig. 6).

Match the number and color indications on the buckles to the shells.

- Adjust the back tension of the Cross Straps.

Ensure that the Cross Buttress is tight against the side of the knee.

Strap chafes can be angled to ensure the straps fit flat against the back side of the leg.

To cut the straps to the desired length, remove the self-fastening hook tab from the end of the strap, cut the strap, and reattach the self-fastening hook tab to the end of the strap (Fig. 7). Take care to not cut any strap too short.

The comfort pads fixed inside the straps may need to be removed for this operation and repositioned after the strap is cut to avoid any interference when the strap is tightened.

For a better fit to the knee, it is possible to bend the upright to conform to the patient's morphology (Fig. 8).

- Place the knee brace on the edge of a table to stabilize the hinge.

- Bend the upright just above or below the hinge.

• Step 3: Check the fit, suspension and pain level prior to adjusting the corrective force

Ask the patient to stand and walk, taking normal steps and looking straight ahead.

Confirm the brace fits and suspends properly on the leg.

Ask the patient if the pain level with the brace on, still in the neutral position, is the same or less than without the brace.

• Step 4: Adjusting the corrective force (Fig. 9)

The Dial on the Cross Strap system is set in a neutral position (with dial pulled out).

- To offload the affected side of the knee, push the dial in, and turn clockwise to increase the corrective force to the desired amount based on the level of pain or support needed.

Start with a half turn of the dial and judge whether the correction level still needs to be adjusted after a few steps.

- Turn the dial counter clockwise to lessen the offload or pull to release all the tension.

• Step 5: Evaluation of pain after adjusting the corrective force

After adjusting the offloading of the brace:

- Ask the patient to walk again and assess his/her pain level.

- Repeat the operation, increasing or reducing the correction until the proper amount of offloading is achieved for the patient.

- If the patient feels any discomfort, reduce the corrective force.

It is recommended to initially wear the brace for only a few hours a day.

It may take several weeks to feel comfortable with the brace on your leg.

Flexion/extension control:

The brace's default extension setting is 0°.

The extension limitations are attached to a plastic tree in the box (Fig. 10).

• Adjusting the extension limitation (Fig. 10):

Extension can be limited to 0°, 5°, 10°, 15°, 20° and 30°.

1. Choose the desired extension limitation on the corresponding plastic piece.

2. Remove the screw located on the hinge side.

3. Flex the hinge slightly and remove the existing extension limitation.

Take care to note the direction this limitation is facing.

4. Insert the desired limitation, numbered end first, with the hook end at the top and facing forward.

Straighten the hinge to full extension to confirm the extension limitation is correctly positioned.

The small hole in each limitation must be aligned and visible through the hole in the hinge, so the screw will thread into the limitation.

5. Reinsert and tighten the screw.

Do a few flexions/extensions to make sure that the limitation is properly locked at the desired angle.

• Adjusting the flexion limitation (Fig. 11):

Flexion can be limited to 0°, 30°, 45°, 60°, 75°, 90° and 110°.

1. Choose the desired flexion limitation on the corresponding metallic piece.

Each limitation has the degree etched into the surface.

2. Remove the two screws from the posterior aspect of hinge cover and take out the default spacer.

3. Insert the limitation with the flat end facing up and position it so that both holes in the limitation are visible through the holes in the cap.

4. Insert and tighten the screws through the cap and into both holes in the limitation.

Do a few flexions/extensions to make sure that the limitation is properly locked at the desired angle.

Warning:

The adjustment of the flexion/extension must be decided and performed by the healthcare professional, not by the patient.

NOTE: It is always better to put your brace on a little too high than too low. Thuasne cannot be held responsible for undesirable effects or injuries resulting from any unsupervised or inadapted adjustments.

How to take off the knee brace:

To take off the knee brace, pull the dial to release tension and unfasten all straps.

Care

Close the self-fastening tabs before washing. Hand wash. Do not use detergents, fabric softeners or aggressive products (products containing chlorine). Squeeze out excess water. Dry away from any direct heat source (radiator, sun, etc.). If the device is exposed to seawater or chlorinated water, make sure to rinse it in clear water and dry it.

• Hinge ②:

The hinge on the brace is pre-lubricated in the factory.

If sand, dirt or water gets inside the hinge, it may require lubricating again.

If you notice the hinge not gliding smoothly, a few drops of a synthetic lubricant can be applied.

Wipe off any excess lubricant before wearing the brace to prevent stains on clothing.

• Straps:

After prolonged use, if the fibers on your strap do not adhere as well to the self-fastening tab, cut the strap shorter so the self-fastening tab adheres to a section of the strap that has fresher fibers.

If this is not possible, you should contact the medical provider who fit your brace.

• Pads:

The brace is lined with padding that provides a comfortable interface between the leg and the shells.

The straps also have pads.

Do not remove these pads from the brace or straps.

Wipe the pads after each use to remove any moisture and let the pads air dry. You can also clean the pads with a mild anti-bacterial soap and rinse them off with fresh water.

Do not wash the pads in a machine and do not dry in a tumble dryer.

Storage

Store at room temperature, preferably in the original packaging.

Disposal

Dispose of in accordance with local regulations.

COMMERCIAL WARRANTY AGREEMENT AND WARRANTY LIMITATIONS

Thuasne offers a free, limited commercial warranty to the user, in the territory where the device was purchased, against defects in manufacturing and workmanship for a period of:

- six months for the straps, padding and condylar pad,

- one year for the shells, buckles and the brace's hinge with its cover.

The limited warranty is effective from the date of purchase of the product by the end-user.

The limited commercial warranty does not apply to any defects in manufacturing and workmanship in case of:

- misuse of the product or any damage occurred by a usage outside the normal and intended use of the product as mentioned in the instructions for use,

- damages occurred as part of attempts to modify the product.

Any deterioration or improper trimming of the product during its modification or adjustment process by the healthcare professional when fitting device is expressly excluded from this warranty.

Any claim for this commercial warranty must be sent by the user to the entity where the product was purchased, which will forward this claim to the corresponding Thuasne entity.

Any warranty claim will first be reviewed by Thuasne to determine if the conditions of the limited warranty are fulfilled and do not fall into one of the cases of exclusion of the commercial warranty.

To benefit from the warranty, the buyer must mandatorily provide an original and dated proof of purchase of the product.

If the conditions of the limited warranty are fulfilled and the claim is made by the user or its legal representative (parents, guardian...) within the warranty delays indicated above, the buyer will get a new substitution product.

It is expressly agreed that this commercial warranty is in addition to the

legal warranties binding the entity which sold the product to the user, in accordance with the applicable local legislation in the country of purchase of the product.

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©Thuasne - 2042101 (2021-05)
**THUASNE DEUTSCHLAND GmbH**
THUASNE UK Ltd
Thuasne SAS - SIREN/RCS Nanterre 542 091 186 - capital 1 950 000 euros