




THUASNE®

Osteoarthritis knee braces prescription guide



If the patient can be fitted with each of these products, please make sure to determine in **ACCORDANCE** with the patient, which product would suit the best their **EXPECTATIONS** and **NEEDS** (Discretion, Desired activities...)

* Max thigh circumference = 71 cm/28"
 (1) OA: Osteoarthritis
 (2) KO: Knee Orthosis / KAFO: Knee-Ankle-Foot orthosis
 (3) Full rigid metallic frame with double upright
 (4) Semi-rigid shells with a metallic upright
 (5) Lamberg, Eric M., Robert Streb, Marc Werner, Jan Kremenc, and James Penna. 2016b. "The 2-and 8-Week Effects of Decompressive Brace Use in People with Medial Compartment Knee Osteoarthritis." *Prosthetics and Orthotics International* 40 (4): 447-453. + Full rigid metallic frame with double upright
 (6) Soft close knitted product without shells (mid stabilization in case of weight). Need to be pulled on the leg upright
 (7) Open product with quick buckles that can be put on the leg without stepping through or pulling on
 (8) Thourme Philippe, et al. "Effect of Unloading Brace Treatment on Pain and Function in Patients with Symptomatic Knee Osteoarthritis: The ROTOR Randomized Clinical Trial." *Scientific Reports* 8 (1): 10519.

Availability of these products might vary from a given country or region to another, as a result of specific local regulatory approval or clearance requirements for sale in such country or region. All the medical devices mentioned in this document are CE marked according to the Regulation 2017/745 on medical devices, they are also UKCA-marked according to S.I.2002/618 regulation on medical devices. Medical devices mentioned in this document are CE class I. Please contact Thuasne® should you need any additional information on devices classification. Please carefully read the instructions for use, indications and contraindications of the products. The full list of medical indications, are available in the instructions for use. Last revision date: 03/2024. Ref.: 2402100. Thuasne SAS - SIREN/RCS Nanterre 542 091 186 - capital 1 950 000 euros - 120, rue Marius AUFAN 92300 Levallois-Perret (France)

Patella Reliever® - Action Reliever® - Rebel Reliever®:
 Thuasne 120, rue Marius AUFAN
 92300 Levallois-Perret
 France

UniReliever™ - SpryStep® KAFO - SpryStep® KO:
 Thuasne Deutschland GmbH
 Im Steinkamp 12
 30938 Burgwedel - Deutschland

Rebel®:
 Thuasne USA
 4615 Shepard Street - Bakersfield, CA 93313
 United-States

MedEnvoy Global BV
 Prinnses Margrietplantsoen 33, Suite 123 - 2595 AM The Hague
 The Netherlands

UK CA All products except Action Reliever® and Patella Reliever®:
 UK Responsible Person (UKRP):
 THUASNE UK Ltd
 Unit 4 Orchard Business Centre
 North Farm Road
 Tunbridge Wells, TN2 3XF,
 United Kingdom