

Contact Information

Clinician Fitter/Assistant/Tech Other: _____

Name: _____

Email: _____ Phone: _____

Ordering Clinician

CPO CO CP Other: _____

Name: _____

Email: _____ Phone: _____

Billing & Shipping

PO#:

Billing Account#: _____

Shipping Account#: _____

Shipping Address: _____

City: _____ State: _____ Zip: _____

Shipping Preference

Ground Next Day AM Next Day PM 2-Day AM 2-Day PM

(If no preference is indicated, this order will be shipped 2 Day PM) Note: We do not ship products directly to patients.

Patient Information

Fit Date: _____ Patient ID: _____

Measuring Information

Measure arm from elbow crease to base of the index finger: S/M up to 14"; L/XL up to 15 1/2".

Low Profile ROM Elbow (Dual Hinges) With Shells And Optional Hand Extension

- Flexion can be limited at 90°, 110° and 120°, and extension can be stopped at 0°, 15°, 30°, 45°, and 90°.
- The screw stops are adjusted using a small flathead screwdriver.
- The hinge can be locked at 90°.

Quantity: _____ S/M _____ L/XL

Install Optional Hand Extension? _____ Yes _____ No

Include Shoulder Sling Strap? * _____ Yes _____ No

Multi-Adjustment ROM Elbow

- Comes with four straps and individual black pads that fit under each strap.
- Tool-less flexion and extension stops are in 10° increments.
- Locking is possible by offsetting the flexion and extension stops at the same degree.

Single Hinge (SH) Model

Quantity: _____ LEFT ARM Length: _____ 15 in _____ 16 in

Quantity: _____ RIGHT ARM Length: _____ 15 in _____ 16 in

Dual Hinge (DH) Model

Length (Indicate quantity): _____ 15 in _____ 16 in

ROM-LBO

- The LBO comes with dual or single uprights, four straps, tool less adjustability, and 8 individual Breeze pads.
- Tool-less flexion and extension stops are in 10° increments.
- The LBO can be locked at 90° for post-surgical protocols.
- Adjustable length 14 1/2" – 17"

Single Hinge (SH) Model

Quantity: _____ LEFT ARM

Quantity: _____ RIGHT ARM

Dual Hinge (DH) Model

Quantity: _____ Universal