

### Contact Information

Clinician  Fitter/Assistant/Tech  Other: \_\_\_\_\_  
 Name: \_\_\_\_\_  
 Email: \_\_\_\_\_ Phone: \_\_\_\_\_

### Ordering Clinician

CPO  CO  CP  Other: \_\_\_\_\_  
 Name: \_\_\_\_\_  
 Email: \_\_\_\_\_ Phone: \_\_\_\_\_

### Billing & Shipping

PO#:

Billing Account#: \_\_\_\_\_  
 Shipping Account#: \_\_\_\_\_

Shipping Address: \_\_\_\_\_  
 City: \_\_\_\_\_ State: \_\_\_\_\_ Zip: \_\_\_\_\_

### Shipping Preference

Ground  Next Day am  Next Day pm  2-Day am  2-Day pm

*(If no preference is indicated, this order will be shipped 2 Day PM) Note: We do not ship products directly to patients.*

### Patient Information

Fit Date: \_\_\_\_\_ Patient ID: \_\_\_\_\_  
 Age \_\_\_\_\_ Weight \_\_\_\_\_ (LBS) Height \_\_\_\_\_ (IN)  
 Male  Female Leg:  Left  Right  
 Diagnosis: \_\_\_\_\_  
 Surgeries (type/date): \_\_\_\_\_

### Brace Configuration

Thigh Shell Length  7"  8"  
 Tibia Shell Length  7"  8"

### Locking Position *(if not indicated, zero degree lock position will be applied)*

0° Locking Position  5° Locking Position

### Size

XS  SM  MD  LG  XL

Sizing	6" above	M-L width	6" below
XS	12 1/2 – 11/2	3" – 3 1/2	11" – 12 1/2
SM	15 1/2 – 18 1/2	3 1/2 – 4"	12 1/4 – 13 3/4
MD	18 1/2 – 21"	4" – 4 1/2	13 1/4 – 15"
LG	21" – 23 1/2	4 1/2 – 5"	14 1/4 – 15 3/4
XL	23 1/2 – 25"	5" – 5 1/2	15" – 17"

### Color

#### Matte Finish

Black (standard)  Grey  White  
 Atlantic (Light Blue)  Red

#### Satin Finish

Lemon  Orange  Lime  
 Fuchsia  Pacific (Dark Blue)


### Accessories

Add Extension Assist Bands/Posts\*  
**Flexion Stop Kit** *(installed upon request)*  
 30°  60°  90°  
 Anti-Migration Silicon Infused Strap Pads\*  
**Undersleeves\***  
 18" Cotton  18" Neoprene

### Measurement Data

<sup>1</sup>These measurements are required to check the accuracy of the patient model submitted, a patient model must be provided for fabrication (cast or scan).

\_\_\_\_\_ Circumference 6 inches above mid-patella  
 \_\_\_\_\_ Medial-Lateral Knee Width (not circumference) at knee center  
 \_\_\_\_\_ Circumference 6 inches below mid-patella



### Special Instructions:

\_\_\_\_\_  
 \_\_\_\_\_

\*Indicates additional charges apply

OF-03 REV E

Received Date *Thuasne USA's shipping department use only*