



THUASNE®

SPINAL BRACING PRODUCTS

*Complete range of
spinal orthoses and
back supports*




**SLEEQ®
MAX**

SLEEQ[®] ADAPT *Premium Step-down LSO*



Individualized support and comfort.



Step-Down System

Modular design easily accommodates progressive therapies.



SLEEQ[®] ADAPT 627

PDAC Approved L0627/
L0642

Low-profile lumbar orthosis
Part # 901700



SLEEQ[®] ADAPT 631

PDAC Approved L0631/
L0648

Lumbar-sacral orthosis
Part # 901800



SLEEQ[®] ADAPT 637

PDAC Approved L0637/
L0650


Lumbar-sacral orthosis with
sagittal-coronal control
Part # 901900

Recommended for: for low back support / sagittal plane motion restriction (and coronal plane motion restriction on 637) / increased intracavity pressure to reduce spine and disc load / acute and/or chronic low back and lower extremity pain (sciatica) / recovery from post-operative spinal procedures or injury.

- Patented InvisAdjust™ universal sizing technology
- Customizable wing insert
- Easy access & modifiable interior support panels
- Angle wings to accommodate waist to hip ratios
- Fits waists 26" to 65" with available EXT extension panel

SLEEQ[®] *The Original Universal LSO*

Patented SLEEQ[®] Invis-Adjust



Universal sizing technology available on all SLEEQ models

Recommended for: Chronic or acute low back pain / degenerative disc disease / bulging or herniated disc / low back sprains or strains / lumbar and trunk muscle weakness / post-op lumbar stabilization / fusion / laminectomy / discectomy / fracture management / prophylactic low back support / decompressive procedures / facet syndrome / spinal stenosis / spondylolisthesis / spondylolysis / spondylosis / radiculopathy / injections.

- Patented InvisAdjust™ universal sizing technology
- Breathable open-cell foam padding
- Self-molding ventilated support panels
- Angle wings to accommodate waist to hip ratios
- Fits waists 26" to 62" with available EXT extension panel



SLEEQ[®] AP

PDAC Approved L0627/
L0642

Low-profile sagittal
(anterior-posterior)
control lumbar orthosis
Part # 901100



SLEEQ[®] AP+

PDAC Approved L0631/
L0648

Enhanced profile sagittal
(anterior-posterior)
control lumbo-sacral
orthosis
Part # 901200



SLEEQ[®] APL

PDAC Approved L0637/L0650

Enhanced profile sagittal-
coronal (anterior-posterior-
lateral) control lumbo-sacral
orthosis
Part # 901300


BOA® SI Belt Ultra low profile SI support

PDAC APPROVED L0621

Recommended for: Sacroiliac joint pain / instability / arthritis / sprains and strains / pelvic girdle and symphysis pubis dysfunction / SI and symphysis pubis issues related to post-partum / fracture management / sciatica.

- Patented Lacer Design Compression System
- Simple pull straps for tightening and adjustment
- Ultra-low-profile design


Part # 51024-x



Size	Waist Circumference
SM	25"-35"
MD	35"-45"
LG	45"-55"

SLEEQ® Accessories

Part #	Description
901070	SLEEQ EXT (extension panel for waist sizes 54"- 62")
901050	SLEEQ ANTERIOR PANEL
901060	SLEEQ LATERAL PANELS (2)
901680	SLEEQ MAX Long Sternal Bar



SLEEQ EXT
Extension Panel

SLEEQ® Hot/Cold Pack

Recommended for: **HEAT:** Relief of muscle spasms, muscle stiffness, joint stiffness, and chronic pain.
COLD: Relief of swelling/inflammation/edema following muscle spasms, sprains, strains, trauma, and acute pain.

- Cold or heat from one pack
- Frost free/drip free
- Skin friendly materials
- Microwave and freezer safe

Part # 1030610020



ATTACHES TO
ANY SLEEQ
BACK BRACE

SLEEQ® Universal TLSOs

SLEEQ® ADAPT 456

PDAC Approved L0456/L0457

Flexible sagittal
(anterior-posterior)
control thoraco-lumbo-sacral
orthosis

Part # 902100

PARTNER  PROGRAM



Recommended for: Stable, osteoporotic vertebral fractures, Kyphosis with poor shoulder posture, Spondylitis, Arthritis, Sagittal plane motion restriction, Post-operative thoracic spinal support, post-injury thoracic spinal support

- Simple Universal Fitting
- Anterior-Posterior Control
- Lightweight & Low Profile
- ComfoGlide™ Tensioning System



- Easy strap adjustment



SLEEQ® MAX

PDAC Approved L0464

Rigid tri-planar control
thoraco-lumbosacral
orthosis

Part # 901600



Over the shoulder or under
axilla strapping options


Recommended for: Moderate to severe sagittal and coronal transverse plane motion restriction / increased intracavitary pressure to reduce spine and disc loading / mechanical or discogenic thoracic lumbar and extremity pain / post-operative spinal stabilization / fusion / degenerative disc disease / post-injury to facilitate healing / spondylolysis / spinal stenosis / spondylolisthesis / spondylosis, laminectomy / discectomy / facet syndrome / spinal muscle weakness / arthritis / low back sprains, strains and prophylactic back support.

Universal Sizing with InvisAdjust™

- Fits 26"-54" waist circumference
- EXT Extension Panel fits up to 62" waist
- Adjustable compression
- Full-Circumference Shell Design
- Modular wings, shells and lateral panels
- Dressing relief cutout

Part # 901600

PARTNER  PROGRAM



Step-Down System

Modular design easily accommodates
progressive therapies.

Universal TLSOs



SLEEQ® FLEX

PDAC Approved L0456/L0457
Flexible sagittal
(anterior-posterior)
control thoraco-lumbo-sacral orthosis
Part # 901400

- Simple Universal Fitting
- Anterior-Posterior Control
- Lightweight & Low Profile
- ComfoGlide™ Tensioning System
- Enhanced Sitting Comfort



Recommended for: Mild to moderate osteoporosis / compression fractures / postural issues and kyphosis / chronic or acute low/mid back pain / osteoarthritis / degenerative disc disease / bulging or herniated disc / low-mid back sprains and strains / lumbar and trunk muscle weakness.

Premium Hyperextension TLSO

Dorso Rigid

PDAC Approved L0472

- Optimized fit and control
- Lightweight telescoping aluminum alloy frame
- Back pad positions horizontally or vertically
- Articulating sternal /pelvic pads

Part # 089001



Size	Pelvic Circumference	Pubic Symphysis to Sternum Height
S	23½" - 29½"	12½"
MS	29½" - 35½"	12½"
M	29½" - 35½"	13¾"
LS	35½" - 41¼"	13¾"
L	35½" - 41¼"	15"
XL	41¼" - 45¼"	16¼"

Recommended for: Compression fractures / traumatic fractures / osteoporotic and metastatic fractures / dorso-lumbar osteoarthritis / immobilization after surgical stabilization of thoracolumbar fractures.



Partner Logo Program

Your logo on every SLEEQ brace*. Ask your Thuasne USA representative for details as well as other product logo opportunities.

